

MU Health Care
Surgeon Picklist Report: Items and Comments

Surgeon: Aggarwal MD, Ajay
Procedure: Arthroplasty Knee Total Replacement

Historical Avg:
Total Cases: 1351

Reporting Class	Item Description	Manuf Number	Open Qty	Hold Qty
<u>Bandages, Dressings, Tape</u>				
	DRSG WND 8X3IN TELFA PLSTR CTTN ABS PD NADH LF STER	1238	1	0
<u>Blades</u>				
	BLADE SAW 3.543X.828IN STRYK PRFRM SER THK.05IN STRNM SGT L S	6221-127-090	1	0
	BLADE SAW 3.937X.984IN STRYK PRFRM SER THK.05IN SGT L S6	6125-127-100	1	0
<u>GLOVES</u>				
	GLOVE SURG 6.5 LF PF SMTH BEAD CUFF STER CRM 11.5MM PROTEXIS	2D72PL65X	1	0
	GLOVE SURG 6.5 LTX PF BEAD CUFF MICRO ROUGHEN N-PYRG STER	30465	1	0
	GLOVE SURG 7 LTX PF SMTH BEAD CUFF UNDERGLOVE STER GRN 288	31270	1	0
	GLOVE SURG 8 LF PF BEAD CUFF STER CRM 12IN PROTEXIS PLISPRN	2D72PL80X	2	0
	GLOVE SURG 8 LTX PF BEAD CUFF TXTR SRFC N-PYRG STER STRW BGL	30480	1	0
<u>IMPLANT ORTHO</u>				
	CEMENT BN SMPX P RADOPQ FD STER	6191-1-001	0	2
<u>LINEN DISPOS</u>				
	COVER 2 TIER BACK HVDYTY LF DISP 72IN TBL STER	DYNJP2TBTC	1	0
	DRAPE SURG 51X47IN ADH STER LF DISP STRDRP CLR	1015	1	0
	DRAPE SURG 72X113IN ANSLP FM HKLP LN HLD R ABS STER LF DISP	DYNJP4015	0	1
	TOWEL 27X17IN CTTN STD PREWASH DELINT ABS BLU DISP SURG STER	MDT2168204	1	0
<u>MED SURG SUPPLY</u>				
	BLADE SAW 3.5X.98IN S6 THK.047IN HVDYTY ORNT TTH SGT L BVL	6125-119-090	0	1
	ELTRD PT RTN 9FT VLAB REM C30- LB PLHSV ACRL FM CORD	E7507	1	0
	MAT ABSORB 40X18IN FLR SGRSF YLW	81860-2	0	0
	MNFLD SUCT NPTN 2 STD 4 PORT WST MGMT SYS NS LF DISP	0702-020-000	1	0
	SET INTPLS SUCT TUBE BN CLN TIP HNDPC	0210-110-000	1	0
	SYRINGE LL 50ML LF STER GRAD N-PYRG DEHP-FR PVC FREE MED	1186000777T	1	0
	SYSTEM BN CMNT REV MX BRKWY NOZ LF	0606-563-000	0	1
	SYSTEM WND IRR SURGIPHOR ANTIMIC STER LF DISP	910100	0	1
	TUBING SUCTION 1/4 X 12 STRAIGHT	OR612	1	0
<u>NEEDLE/SYRINGE</u>				
	SYRINGE 60 ML LF STER IRR FNGR FLNG PLBL BLB DISP	DYND20125	2	0
<u>ORTH SOFTG/SUP</u>				
	CRADLE XTRMT PED LF 24X5X3IN ARM LMI	ALP 303	0	0
<u>PACKS MED SUP</u>				
	KIT CUSTOM ROOM TURNOVER OR	DYKC1547	0	0
	PACK SURGICAL T/K STER DISP	DYNJ84110C	1	0
	TOTAL KNEE ACCESSORY PACK	DYNJ45987M	1	0
<u>PEEL PACK</u>				
	BATTERY SYSTEM 9 LARGE (STRYKER) MOI		3	0
	PP ELEVATOR FREER MOI		1	0
	WS BASIN MOI		1	0
	WS BUMP LARGE MOI		1	0
	WS LIGHT HANDLE		1	0
<u>TRAYS, MEDICAL-SURGICAL</u>				
	TRAY CEMENT GUN REVOLUTION STRYKER		0	1
	TRAY MOI TOTAL KNEE		1	0
	TRAY TOTAL JOINT SYSTEM 9 STRYKER- MOI		1	0
	WS STULBERG LEG HOLDER MOI		1	0
<u>WOUND CLOSURE</u>				
	ADH SKNCLS 2-OCTYL CYNCRLT DRMBND PRINEO LQ 22CM	EF70401	0	1
	SUTURE PDO STRATAFIX 1 CTX 36CM ABS KNTLS TISS CNTRL 2 ARM	SXPP2B405	1	0
	SUTURE PLGCPRN STRATAFIX MNCRL + 3 PS SPRL 12IN UD ABS	SXMP2B412	1	0
	SUTURE STRATAFIX PDS + 2 CT1 SPRL 9IN ABS	SXPP2B408	1	0
	SUTURE VICRYL 2-0 CT1 18IN CR BRD 8 STRN COAT UD ABS	J839D	1	0
	SUTURE VICRYL+ 1 CTX 18IN CR BRD 8 STRN ANBCTRL VIOL ABS	VCP765D	1	0

Surgeon: Aggarwal MD, Ajay
Procedure: Arthroplasty Knee Total Replacement

Historical Avg:
Total Cases: 1351

Procedure Specific Preference Card Comments:

PROCEDURE NAME: Arthroplasty Knee Total Replacement (NEW VERSION)
DR. AGGARWAL
GLOVE: 8 Biogel X2 & 8 White
Katie: 7 Green, 6.5 Biogel, 6.5 White

ROOM SETUP:
Skytron bed w/ arm boards
Bovie, tourniquet, & suction always on left side
stryker helmets & batteries for helmets

ESU:
60/60

POSITIONING:
supine, tourniquet, Stulberg knee positioner
may need blanket bump

PREP/CLIP/FOLEY:
tourniquet, 1015
sterillum, chloroprep x2, then 1 on knee,
prep tourniquet to toes

DRAPING:
down sheet, sm. split sheet, stockinette, coban, extremity sheet, 6650 drape

SPECIMEN:

SPECIAL MEDICATION INSTRUCTIONS:
sodium chloride 3L w/ 3 grams ancef, if allergic use plain saline
Call pharm @ 43055 for Ortho Cocktail
1L NACL

X-RAY:
None usually

IMPLANTS:
Stryker implants

SPECIFIC SURGEON PREFERENCE:
for infected case- zimmer palacoec R cement
Simplex with tobra

Dressing:
Exofin strip, telfa, 4x4's, ABD Pads, webril, Ace Wrap, knee immobilizer

Other General Provider Comments:

MU Health Care
Surgeon Picklist Report: Items and Comments

Surgeon: Aggarwal MD, Ajay
Procedure: Arthroplasty Hip Total Replacement

Historical Avg:
Total Cases: 1552

Reporting Class	Item Description	Manuf Number	Open Qty	Hold Qty
<u>Bandages, Dressings, Tape</u>				
	SPONGE GZE 4X4IN CTTN 12 PLY WVN FOLD EDG LF STER	PRM21426	1	0
	WRAP COMPR 5YDX4IN CBN POR SLFADH ELC LTWT HND TEAR LF STER	2084S	2	0
<u>Blades</u>				
	BLADE SAW 77.6X11.2X.77MM LNG OFST SS RECIPROCATATE HVYDTY	277-96-325	1	0
<u>GLOVES</u>				
	GLOVE SURG 6.5 LF PF SMTH BEAD CUFF STER CRM 11.5MM PROTEXIS	2D72PL65X	1	0
	GLOVE SURG 6.5 LTX PF BEAD CUFF MICRO ROUGHEN N-PYRG STER	30465	1	0
	GLOVE SURG 7 LTX PF SMTH BEAD CUFF UNDERGLOVE STER GRN 288	31270	1	0
	GLOVE SURG 8 LF PF BEAD CUFF STER CRM 12IN PROTEXIS PLISPRN	2D72PL80X	2	0
	GLOVE SURG 8 LTX PF BEAD CUFF TXTR SRFC N-PYRG STER STRW BGL	30480	1	0
<u>LINEN DISPOS</u>				
	COVER 2 TIER BACK HVYDTY LF DISP 72IN TBL STER	DYNJP2TBTC	1	0
	DRAPE SURG 51X47IN ADH STER LF DISP STRDRP CLR	1015	1	0
	DRAPE SURG 72X113IN ANSLP FM HKLP LN HLDL ABS STER LF DISP	DYNJP4015	1	0
<u>MED SURG SUPPLY</u>				
	DRAPE SURG 43528 77X53IN FNFLD ABS REINF STER LF DISP PRXM	DYNJP2416	0	1
	ELTRD PT RTN 9FT VLAB REM C30- LB PLHSV ACRL FM CORD	E7507	1	0
	HOOD SURGEON FLYTE STERI-SHIELD SURGICOOL FACE SHIELD VIRAL	0408-800-000	4	1
	MAT ABSORB 40X18IN FLR SGRSF YLW	81860-2	0	0
	MNFLD SUCT NPTN 2 STD 4 PORT WST MGMT SYS NS LF DISP	0702-020-000	1	0
	SET INTPLS SUCT TUBE BN CLN TIP HNDPC	0210-110-000	1	0
	SYRINGE LL 50ML LF STER GRAD N-PYRG DEHP-FR PVC FREE MED	1186000777T	1	0
	SYSTEM WND IRR SURGIPHOR ANTIMIC STER LF DISP	910100	0	1
<u>NEEDLE/SYRINGE</u>				
	SYRINGE 60 ML LF STER IRR FNGR FLNG PLBL BLB DISP	DYND20125	1	0
<u>ORTH SOFTG/SUP</u>				
	CRADLE XTRMT PED LF 24X5X3IN ARM LMI	ALP 303	0	0
	PILLOW LG CNTR 25X18X6IN FM ABD LF DISP	APG201	1	0
<u>PACKS MED SUP</u>				
	KIT CUSTOM ROOM TURNOVER OR	DYKC1547	0	0
	PACK SURG CMP TOTAL HIP CUST LTX	DYNJ84336B	1	0
	PACK SURG T/H STER DISP	DYNJ45991N	1	0
<u>PEEL PACK</u>				
	BATTERY SYSTEM 9 LARGE (STRYKER) MOI		2	0
	WS BASIN MOI		1	0
	WS BUMP LARGE MOI		1	0
	WS LIGHT HANDLE		1	0
<u>TRAYS, MEDICAL-SURGICAL</u>				
	TRAY MOI CURETTE		1	0
	TRAY MOI TOTAL HIP		1	0
	TRAY RETRACTORS TOTAL HIP MOI		1	0
	TRAY TOTAL JOINT SYSTEM 9 STRYKER- MOI		1	0
<u>WOUND CLOSURE</u>				
	ADH SKNCLS 2-OCTYL CYNCRLT DRMBND PRINEO LQ 22CM	EF70401	1	0
	SUTURE PDO STRATAFIX 1 CTX 36CM ABS KNTLS TISS CNTRL 2 ARM	SXPP2B405	1	1
	SUTURE PLGCPRN STRATAFIX MNCRL + 3 PS SPRL 12IN UD ABS	SXMP2B412	1	0
	SUTURE STRATAFIX PDS + 2 CT1 SPRL 9IN ABS	SXPP2B408	1	0
	SUTURE VICRYL+ 2-0 CT1 18IN CR BRD 8 STRN ANBCTRL UD ABS	VCP839D	1	1

Surgeon: Aggarwal MD, Ajay
Procedure: Arthroplasty Hip Total Replacement

Historical Avg:
Total Cases: 1552

Procedure Specific Preference Card Comments:

PROCEDURE NAME: Arthroplasty Hip Total Replacement (NEW VERSION)
DR. AGGARWAL
GLOVE: 8 Biogel X2 & 8 White
Katie: 7 Green, 6.5 Biogel, 6.5 White

ROOM SETUP:

Bovie & Suction always on left side of pt.
1 padded mayo stand w/ foot pedal to change height
Sequential device for non-operative leg & operative/ non- operative
ted hose up to knees
peg board with 4 posts, gel ax roll

ESU:
60/60

POSITIONING:

skytron bed, pegboard w/ gelpad. Secure peg board w/ peg board clamps
(1 @ head, 1 on the other side @ bottom)
Axillary roll
Pillows between arms
Smaller mayo stand covered w/ pillow, secure pillow w/ tape
Make sure edges are padded
Scd on unaffected side

PREP/CLIP/FOLEY:

foley-yes
1010 and 1015 to drape before prep
sterillium wipe before prep to ankle
chloraprep x3 and x1 chloraprep

DRAPING:

down sheet, up sheet, sticky u x2, towel, stockinette, hip drape, cloth stockinette coban x2,
marking pen, ioban x2,
Irrigation pouch is clamped (glove change) to the drape to put the bovie & power irrigator in.

SPECIAL MEDICATIONS:

Cefazolin (Ancef)- 3gms in 3000cc NACL, if allergic use plain saline
NACL on field
Call pharm @ 43055 for Ortho Cocktail

SPECIMEN:

Femoral head for Gross only

SCRUB:

You will need 3 mayo stands- Pad one w/ pillow & coban, drape sterile w/ mayo cover & towels
use small strips of ioban to secure towels to the mayo (He uses this to prop up the operative
leg)
Use one #5 ethibond cuts off needle & tags for the acetabulum he starts reaming 6 under what
is templated
(putting in a 56, will start w/ 50)
Dr. will stop reaming 1 under what he will be inserting
(putting in a 56 will ream to 55)

Dressing:

Exofin strip, telfa, 4x4's, abd, medipore tape, abduction pillow

long debakey forceps & regular weitlaners

have sofieled retractor, hibbs retractors, dull weitlanders, in room (all in center core)

For hip resurfacing will use:

Surgeon: Aggarwal MD, Ajay
Procedure: Arthroplasty Hip Total Replacement

Historical Avg:
Total Cases: 1552

BHR- 60 cc cath-tip syringe
-hibbs with a hook
-2 narrow hohmann retractors
-system 6 collett set (stryker)
-power irrigator (ref 210-110 stryker) interpulse handpiece
-1 box simplex with tobramycin
-zimmer mixing bowl
-3 marking pens

SMITH AND NEPHEW TOTAL HIP TRAYS:

Anthology tray 1, Anthology tray 2, R3 main instrument, R3 reamers, R3 trial shells, MIS tray
(need 6 trays)

Other General Provider Comments:

MU Health Care
Surgeon Picklist Report: Items and Comments

Surgeon: Aggarwal MD, Ajay
Procedure: Arthroscopy Hip

Historical Avg:
Total Cases: 619

Reporting Class

Item Description	Manuf Number	Open Qty	Hold Qty
<u>Bandages, Dressings, Tape</u>			
DRSG TRNS 12X8IN PU ADH HPOAL WTPRF TGDRM STER LF	1629	2	0
SPONGE GZE 4X4IN CTTN 12 PLY WVN FOLD EDG LF STER	PRM21426	1	0
WRAP COMPR 5YDX4IN CBN POR SLFADH ELC LTWT HND TEAR LF STER	2084S	2	2
<u>Blades</u>			
BLADE SHAVER 4MM FRML LNG ANG TOMCAT STER	0385-545-100	1	0
<u>GLOVES</u>			
GLOVE SURG 6.5 LF PF SMTH BEAD CUFF STER CRM 11.5MM PROTEXIS	2D72PL65X	1	0
GLOVE SURG 6.5 LTX PF BEAD CUFF MICRO ROUGHEN N-PYRG STER	30465	1	0
GLOVE SURG 7 LTX PF SMTH BEAD CUFF UNDERGLOVE STER GRN 288	31270	1	0
GLOVE SURG 8 LF PF BEAD CUFF STER CRM 12IN PROTEXIS PLISPRN	2D72PL80X	2	0
GLOVE SURG 8 LTX PF BEAD CUFF TXTR SRFC N-PYRG STER STRW BGL	30480	1	0
<u>IMPLANT ORTHO</u>			
ANCHOR SUT 1.4MM 25D NANOTACK FLEX INSRTR 2 STG FIX PILT H	CAT01858	1	2
DEVICE SUT 1.5MM CRSNT MNGR REACH SM FTPrNT NANOPASS ACETAB	CAT02298	1	0
<u>LINEN DISPOS</u>			
DRAPE SURG 125X83IN 2 INCISE FILM PCH ISL ADH STRP STER STR	6617	1	0
DRAPE SURG 23X17IN 2 INCISE FILM ANTIMIC STER IOBN	6650EZ	1	0
<u>MED SURG SUPPLY</u>			
APL 70% ISPRP 2% CHG 26ML 8.4X8.4IN CHLRPRP HI-LT ORNG SCRB	930815	2	0
BAG 30X30IN BAND EQP	01-3030	1	0
BIT DRILL NANOTACK FLEX ANATM LABRUM RSTR SYS	CAT01859	1	0
BURR SHAVER 5.5MM LNG 8 FLUT RND F SER	0385-950-008	1	0
CANNULA ARTHRO 5MM 165MM FLOWPORT II STRYK PLMR HIP OBT	CAT02438	1	0
CANNULA ORTHO 7 8 9 8MM TRNSPT	CAT00224	1	0
DEVICE CLOSURE CAPSULE SLINGSHOT 45 DEGREE	CAT02589	1	0
ELECTRODE ESURG 110D DOME 3.4MM ARCTIC VAPR HIP SUCT LOPRO	225029	1	0
KIT ARTHRO 17GA HIP NITINOL PRTL ENTRY 2 NDL SYR GWIRE 30CC	CAT00241	1	0
LINER BT UNIV DISP	5937DZ	1	0
MAT ABSORB 40X18IN FLR SGRSF YLW	81860-2	0	0
MNFLD SUCT NPTN 2 STD 4 PORT WST MGMT SYS NS LF DISP	0702-020-000	1	0
PAD PSTN 11IN 9IN 1.5IN ACT HL HIP PERINEUM SUPN POST DISP	72200634	1	0
SUTURE ZIPLINE 36IN COBRAID B/W NABSB	CAT01586	2	2
TUBING SUCTION 1/4 X 12 STRAIGHT	OR612	1	0
<u>PACKS MED SUP</u>			
PACK SURG HIP ARTHRO STER DISP	DYNJ81884B	1	0
<u>PEEL PACK</u>			
BATTERY SYSTEM 9 LARGE (STRYKER) MOI		1	0
PP STRYKER HIP BLADE SHARP TIP CVD MOI		0	1
WS SHAVER ARTHROSCOPY STRYKER MOI		2	0
<u>TRAYS, MEDICAL-SURGICAL</u>			
PIVOT HIP STRYKER TRAY MOI		1	0
TRAY ARTHROSCOPY HIP REPAIR EXTRAS MOI		0	1
TRAY CAMERA ARTHREX W/COUPLER MOI		1	0
TRAY HIP ARTHROSCOPY REPAIR MOI		1	0
TRAY SAW SABO SYSTEM 9 STRYKER- MOI		1	0
TRAY SCOPE HIP SINGLE 4MM 70DEG STRYKER MOI		1	0
TRAY STRYKER ARTHROSCOPY HIP		1	0
<u>WOUND CLOSURE</u>			
STAPLER SKN 3.9MM 360D REG RECT 35 CNT ROT HEAD RCHT STER LF	PRR35	1	0
SUTURE ETHILON 3-0 FS1 18IN MFL BLK NABSB	663G	1	0

Surgeon: Aggarwal MD, Ajay
Procedure: Arthroscopy Hip

Historical Avg:
Total Cases: 619

Procedure Specific Preference Card Comments:

PROCEDURE NAME: Arthroscopy Hip
DR. AGGARWAL
GLOVE: 8 Biogel X2 & 8 White
Katie: 7 Green, 6.5 Biogel, 6.5 White

ROOM SETUP:

Hana table large white foam perineal pad
Depuy Vapor w/ foot pedal, arthrex pump, stryker shaver box x2 with foot pedals
Tower, fluid pump @ 50
suction on non operative side
C-arm between legs & endocam #3 flouro on middle monitor (DUAL MONITORS ANGLED BY NON OP)
Printer set at 6/2; will always record, burn to dvd
C-arm monitor to foot of bed

ESU:

No (Uses Depuy Vapor at 60/60)

POSITIONING:

Hana Table w/ the white post, coban to feet and over boot, slide pt down to fit snugly
against padded perineal post
Operative side arm across chest w/ blanket and pink foam pad, secured w/ coban, non-operative
side arm to arm board, bair hugger

PREP/CLIP/FOLEY:

No pre-drape
Sterillium pre prep
Chloraprep lower chest to knee x3

DRAPING:

towels and stapler, Ioban shower curtain drape, Lap Chole drape, 2 mayo stands - 1 will be
place on the non-op side over the pt. place all of the cords & tubes camera etc 1 suction
connected to bag

SPECIMEN:

None

SPECIAL MEDICATION INSTRUCTIONS:

3cc of 1:1000 epi/ 3L bag of LR , Start W/ 6 Bags
30cc of 0.25% Bupivacaine with epi at end of case

SPECIFIC SURGEON PREFERENCE:

all monitors to opposite side
circ let record at start of case

If doing proximal femoral osteotomy need:

Synthes compact battery drill
Synthes bone forceps
Arthrex bio suture inst- hold

Note traction time on the OR record & operative time information form
record whole process *burn dvd and label w pt label
Make 2 copies of pics

Starts w/ 70 degree lens, gray cannulas x2, clear cannula x1, uses 2 spinal needles & 2
guidewires & the bridge that goes on the lens.

Dressing:

adaptic, 4x4's, tegaderm

Other General Provider Comments: